

Ministat® 230

Ministat® 230 cycling a 500 mL DWK flat bottomed bottle reactor

Requirement

This Case Study demonstrates the temperature control capabilities of the process temperature when a Ministat 230 is connected with a 500 mL DWK flat bottomed bottle reactor.

Method

The 500-mL flat bottomed bottle reactor was connected to Ministat® 230 using Tygon A-60-G hoses. The thermofluid used in the system was "M40.165/220.10", reactor was filled with 400 mL Isopropyl Alcohol. "Process" control was carried out. Magnetic stirrer speed was set to 200 rpm.

Setup details

Temperature range: -40°C...+200°C
 Cooling power: 0.42 kW @ +20°C
 0.38 kW @ 0°C
 0.25 kW @ -20°C
 Heating power: 2.0 kW
 Hoses: Tygon A-60-G
 HTF: M40.165/220.10
 Reactor: DWK 500 mL flat bottomed bottle reactor
 Reactor content: 400 mL Isopropyl Alcohol
 Stirrer speed: 200 rpm
 Control: process



Results

Performance:

The graphic shows the control of this simple bench-top application.

Start T	End T	Time Taken	Av. Ramp Rate
+20°C	50°C	31 minutes	1.0 K/min
+50°C	-29°C	T-min	xxxx
-29°C	+20°C	36 minutes	1.4 K/min

