

Unistat® 930w

Controlling a Buchi Glas Uster CR101, a 100-litre GLSS reactor from 20 °C to -40 °C

Requirement

The graphic demonstrates the capability of a Unistat 930w to cool a 100-litre reactor from 20 °C to -40 °C.

Method

The Unistat and reactor are connected using two 2-metre insulated metal hoses. The reactor is filled with 75 litre of "M90.055.03", a Huber supplied silicon based HTF.

Results

The process temperature reaches its target of -40 °C within 43 minutes representing a ramp-rate > 1.4 K/min.

Setup details

Unistat® 930w & 100-litre Buchi Glas Uster «chemReactor» CR101

- Temperature range: -90...200 °C
- Cooling power: 20 kW @ 0...-40 °C
15 kW @ -60 °C
- Heating power: 24 kW
- Hoses: 2x1.5 m; M38x1.5 (#6656)
- HTF: DW-Therm (#6479)
- Reactor: 100-litre glass-lined (enameled) steel reactor
- Reactor content: 75 litre M90.055.03 (#6259)
- Stirrer speed: 80 rpm
- Control: process

