

Unistat® 910w

Cooling a Radleys 10-litre reactor from 180 °C to 20 °C

Requirement

This graphic shows the performance of a Unistat 910w cooling a Radleys 10-litre glass reactor from 180 °C to 20 °C under process control.

Method

The Unistat and reactor are connected using two 1.5-metre insulated metal hoses. The reactor is filled with 7.5 litre of "M90.055.03", a Huber supplied silicon based HTF.

Results

The „internal“ (jacket) temperature quickly cools to -28 °C pulling the process temperature to the set-point.

The process temperature ramps through 160 K (180 °C to 20 °C) in 32 minutes, a process ramp rate of 5 K/min.

Setup details

Unistat® 910w & Radleys reactor

Temperature range: -90...250 °C
 Cooling power: 5.2 kW @ 250...-20 °C
 4.7 kW @ -40 °C
 Heating power: 6.0 kW
 Hoses: 2x1.5 m; M30x1.5 (#6386)
 HTF: DW-Therm (#6479)
 Reactor: 10-litre glass reactor (#6259)
 Reactor content: 7.5 litre M90.055.03
 Stirrer speed: 200 rpm
 Control: process

